

Peri-Strips Dry®

With Veritas® Collagen Matrix • Staple Line Reinforcement
Avec matrice de collagène Veritas® • Renforcement de ligne d'agrafes
Mit Veritas® Collagen-Matrix • Klammernaht-Verstärkung
Con matrice di collagene Veritas® • Rinforzo della linea di sutura
Con matriz de colágeno Veritas® • Refuerzo de la línea de grapado
Med Veritas® Collagen Matrix • Nietertlijnversterking
Med Veritas® Kollagen -matrix • Hæftesømsforstærkning
Med Veritas® Kollagenmatrix • Agrafflinjeförstärkning
Med Veritas® kollagen-matrise • stiftelinjeforstærkning
Veritas® Kollagen Matriksli • Zimba Hattu Kuvvetlendirici
με μήτρα κολλαγόνου Veritas® • Υλικό ενίσχυσης γραμμής συνδέτηρων συρραπτικού
Veritas® kollageenimatriksilla • Hakasrivin vahvike
Com Matriz de Colagénio Veritas® • Reforço da Linha de Agrafa
Cu matrice de colagen Veritas® • Intărirea liniei de sutură
S kollagenovou maticí Veritas® • Zpevnění linie svorek
Veritas® kollagénmátrixszal • Kapocsosor-megerősítő
Z kolagenową macierzą Veritas® • Wzmocnienie linii szwu mechanicznego

INSTRUCTIONS FOR USE	4
MODE D'EMPLOI	7
GEBRAUCHSANLEITUNG	10
ISTRUZIONI PER L'USO	14
INSTRUCCIONES DE USO	17
GEBRUIKSAANWIJZING	20
BRUGSVEJLEDNING	23
BRUKSANVISNING	26
BRUKSANVISNING	29
KULLANMA TALIMATI	32
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	35
KÄYTTÖOHJEET	39
INSTRUÇÕES DE UTILIZAÇÃO	42
ISTRUCȚIUNI DE FOLOSIRE	45
NÁVOD K POUŽITÍ	48
HASZNÁLATI ÚTMUTATÓ	51
INSTRUKCJA UŻYCIA	54

Peri-Strips Dry® with Veritas® Collagen Matrix Components

Peri-Strips Dry® with Veritas® Collagen Matrix Figures 1-7


FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4


FIGURE 5


FIGURE 6


FIGURE 7

A

Loading Unit
Unité de chargement
Gewebebausatz
Unità di caricamento
Unidad de carga
Laadeinheit
Klargöringsenhet
Laddningsenhet
Ladeinheit
Yükleme Ünitesi
Μονάδα φόρτωσης
Lataaja
Unidade de Carga
Unitate de încărcare
Nakládací jednotka
Betöltőegység
Wklad

B

Pericardium Buttress Strips
Bandes de renfort de péricarde
Pericardium-Streifen
Strisce di materiale di sostegno pericardico
Tiras de soporte del pericardio
Hartzakje Steunstrips
Peri-strips (strimler til perikardium) m. forstærkning
Stödremсор av perikardium
Perikardium forstöttningsstrips
Perikard Destek Stripleri
Ταινίες υλικού ενίσχυσης από περικάρδιο
Vahvikkeet sydänpuussia varten
Tiras de Apoio de Pericárdio
Benzi de susținere din pericard
Zpevňovací proužky perikardu
Perikardium támszalagok
Paski wzmacniające z osierdzia

C

Plastic Sheath
Gaine en plastique
Plastikhülle
Guaina in plastica

Funda de plástico
Plastic huls
Plastichylster
Plasthölje
Plasthylse
Plastik Kilf
Πλαστικό θηκάρι
Muovikotelo
Revestimento de Plástico
Folie de plastic
Plastový plášť
Műanyag tok
Plastikowa koszulka

D

Assembly Clip
(ENDO GIA UNIVERSAL Only)
Clip d'assemblage
(ENDO GIA UNIVERSAL seulement)
Bausatzklammer
(nur ENDO GIA UNIVERSAL)
Fermo del gruppo
(solo ENDO GIA UNIVERSAL)
Clip de montaj
(sólo ENDO GIA UNIVERSAL)
Assemblageklem
(Alleen ENDO GIA UNIVERSAL)
Sætklemme
(Kun ENDO GIA UNIVERSALKLEMMEN)
Monteringsklämma
(Endast ENDO GIA UNIVERSAL)
Sammensettingsklips
(ENDO GIA bare UNIVERSAL)
Kurma Klipsi
(Sadece ENDO GIA EVRENSEL)
Κλιπ διάταξης
(ENDO GIA UNIVERSAL μόνο)
Asennusliitin
(Koskee vain ENDO GIA UNIVERSAL -mallia)
Clip de Montagem
(Somente ENDO GIA UNIVERSAL)
Clemă de asamblare
(Numai pentru tipul ENDO GIA UNIVERSAL)

Svërka souprav (pouze ENDO GIA UNIVERSAL)
Össetartó kapocs (csak az ENDO GIA UNIVERSAL esetében)
Zacisk montazowy (tylko dla modeli ENDO GIA UNIVERSAL)

E

Peri-Strips Dry Gel
Gel Peri-Strips Dry
Peri-Strips Trockengel
Gel Peri-Strips Dry
Gel Peri-Strips Dry
Peri-Strips Dry Gel
Peri-Strips Dry gel
Peri-Strips torr gel
Peri-Strips Dry Jel
Γέλη Peri-Strips Dry
Peri-Strips kuivageeli
Gel Peri-Strips Dry
Gel Peri-Strips Dry
Gel pro proužky Peri-Strips Dry
Peri-Strips Dry gél
Zel Peri-Strips Dry

F

Gel Tip
Extrémité du gel
Gelspitze
Punta con gel
Punta del tubo de gel
GelTip
Gelspids
Gelspets
Gelspiss
Jel Űcu
Ρύγχος γέλης
Geelikärki
Ponta de Gel
Aplicator gel
Špička pro aplikaci gelu
Gélhegy
Końcówka zelowa

SYMBOL DEFINITIONS:

20°C/68°F  25°C/77°F

Store at a controlled room temperature.

Keep away from heat. Do not use if heat indicator is red.

Do not reuse

Consult *Instructions for Use*

STERILE

Sterile

STERILE EO

Sterilized using ethylene oxide

Do not use if the product sterilization barrier or its packaging is compromised.

STERILE

Sterilized using steam heat

NaOH

This product is treated with sodium hydroxide.

BOVINE

This product is derived from USDA inspected cattle.

MADE IN THE U.S.A.

Made in the USA

Rx Only

CAUTION: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

REF

Catalog number

Use by date

LOT

Lot number

Compatible stapler models

REMOVE BEFORE USE

Remove before use

Manufacturer

EC REP

Authorized Representative in the European Community

PN

SSI part number

TN

SSI tracking number

IC

SSI internal code

ID

SSI internal code

DESCRIPTION:

Peri-Strips Dry® with Veritas® Collagen Matrix (PSD-V) is prepared from dehydrated bovine pericardium procured from cattle under 30 months of age in the United States.

One (1) tube of PSD Gel (Gel) is provided for every two (2) pouches of PSD-V. The Gel is used to create a temporary bond between the PSD-V buttress and the surgical stapler jaws until the stapler is positioned and fired. Each PSD-V loading unit and each Gel tube is packaged sterile in a separate pouch.

INDICATIONS FOR USE:

PSD-V is intended for use as a prosthesis for the surgical repair of soft tissue deficiencies using surgical staplers when staple line reinforcement is needed.

PSD-V can be used for reinforcement of staple lines during lung and bronchus resections and during bariatric surgical procedures.

PSD-V can be used for reinforcement of staple lines during gastric, small bowel, mesentery, colon, and colorectal procedures.

CONTRAINDICATIONS:

The use of PSD-V is contraindicated in patients with known sensitivity to bovine material.

ADVERSE REACTIONS:

As with any surgical procedure, adverse reactions are possible and include but are not limited to: infection, rejection, erosion, and allergic reaction.

WARNINGS:

Do not re-sterilize. Resterilization may cause changes to the tissue and negatively impact functionality of the device.

Discard all open unused components.

Do not use product if there is damage to the pouch or seals.

Ensure the staple line is completely covered with the buttress, or inadequate coverage after firing may result.

PSD-V is not designed, sold, or intended for use except as indicated; doing so may result in surgical complications.

Synovis products differ; substitution of one product for another product may reduce product performance.

Do not use PSD-V if the product has been exposed (1) to solutions above room temperature, (2) to chemicals, antibiotics, or other substances other than specifically addressed in the instructions, (3) or if heat indicator has been activated.

CAUTIONS:

Use care when removing loading unit components from the stapler to prevent buttress dislodgement.

Do not get the buttress wet before applying Gel, or the buttress may not adhere to the stapler properly.

Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the buttress may not adhere to the stapler properly.

Final tissue compression, including PSD-V, must meet the range specified by the stapler manufacturer; this is especially important if staple firings are overlapped. PSD-V increases the total thickness of the area stapled by 0.4 mm - 1.2 mm (0.016" - 0.048").

Follow *Instructions for Use* supplied by the stapler manufacturer. Do not use PSD-V contrary to the stapler manufacturer's instructions.

The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.

INSTRUCTIONS FOR USE:

LOADING STAPLER FOR ALL STAPLER MODELS

Note: *Each model of PSD-V has been designed specifically for the stapler models indicated on the label; verify that the correct model of PSD-V has been selected.*

Note: *Loading technique for PSD-V varies. Follow the appropriate technique as indicated below.*

1. Open the outer PSD-V pouch.
2. Aseptically remove the inner pouch and place in the sterile field.
3. Inspect the pouch. Do not use if the pouch is damaged or if the seals are not intact.

4. Open the inner PSD-V pouch and remove the PSD-V loading unit by using atraumatic techniques.
5. Inspect the Gel pouch. Do not use if the pouch is damaged or if the seals are not intact.
6. Open the Gel pouch and aseptically place Gel in the sterile field.
7. Break off the metal tip of the Gel tube.
8. Attach the plastic Gel Tip to the Gel tube by pushing it down over the top of the Gel tube.
9. Apply a continuous thin Gel bead onto each buttress strip (See Figure 1).

Caution: Do not get the buttress wet before applying Gel or the buttress may not adhere to the stapler properly.

10. Ensure that the jaws of the stapler are clean and dry before inserting the PSD-V loading unit.

Caution: Final tissue compression, including PSD-V, must meet the range specified by the stapler manufacturer. This is especially important if stapler firings are overlapped. PSD-V increases the total thickness of the area stapled by 0.4mm - 1.2mm (0.016" - 0.048").

11. Identify the anvil (ANV) and cartridge (CART) sides of the PSD-V loading unit.

Caution: The cartridge and anvil loading units differ; substitution of one unit for another may interfere with closing and firing of the stapler.

Note: *For ENDO GIA™ UNIVERSAL stapler models skip to step 19.*

LOADING STAPLER FOR ALL STAPLER MODELS EXCLUDING ENDO GIA™ UNIVERSAL MODELS

12. Position the open stapler onto the stapler loading unit (See Figure 2).
13. Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the buttress may not adhere to the stapler properly.

Note: *The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.*

14. Close the stapler.
15. Remove the plastic sheath leaving the foam spacer between the stapler jaws (See Figure 3).
16. Discard plastic sheath and remove any excess Gel from stapler.

Note: PSD-V can be used immediately or allowed to remain between the stapler jaws until surgery.

Note: For optimal results the prepared PSD-V/stapler should be used within 60 minutes.

17. Just prior to using the stapler, open the stapler and remove the foam spacer. (See Figure 4).
18. Visually inspect each buttress strip to ensure they are on the stapler jaws and cover the staple line.

Note: Skip to step 27.

LOADING STAPLER FOR ENDO GIA™ UNIVERSAL MODELS

19. Position the open stapler onto the stapler loading unit (See Figure 2).
20. Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the strips may not adhere to the stapler properly.

Note: The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.

21. Close the stapler.
22. Apply light pressure to both sides of the stapler jaws to hold stapler into closed position.
23. Slide the assembly clip over the stapler jaws (See Figure 5).
24. Ensure that both stapler jaws are held firmly in place by the assembly clip and that the yellow guard, if present, is positioned between the assembly clip on the cartridge side.

Note: PSD-V can be used immediately or allowed to remain between the stapler jaws until surgery.

Note: For optimal results the prepared PSD-V/stapler should be used within 60 minutes.

25. Immediately before using, remove stapler loading unit (See Figure 6).
26. Visually inspect each buttress strip to ensure they are on the stapler jaws and cover the staple line.

IMPLANTING PSD-V INSTRUCTIONS FOR ALL STAPLER MODELS

27. Follow *Instructions For Use* supplied by the stapler manufacturer.
28. If necessary, cut the end of the PSD-V to remove the dissected tissue (See Figure 7).
29. Discard any open PSD-V and Gel pouches. These cannot be re-sterilized or reused.

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written, or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device, which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Peri-Strips Dry is a registered trademark of Synovis Life Technologies, Inc., (US Patent #5,752,965; 6,312,474; 6,652,594; Additional Patents Pending)
A Division of Synovis Life Technologies, Inc.