
eng

Baxter

REF ADS201865

EASYGRIP FLO-41 Precision MIS Delivery System

INSTRUCTIONS FOR USE

Device Description:

The EASYGRIP FLO-41 Precision MIS Delivery System ("EASYGRIP FLO-41 System")

is a sterile, single-use device that consists of two components:

- 1 x applicator device with a 41cm long cannula (5mm outer diameter)
- 1 x empty 1.5 mL syringe.

The EASYGRIP FLO-41 System is compatible with FLOSEAL Hemostatic Matrix (hemostatic agent).

Sterile and non-pyrogenic. Sterilized using irradiation.

Not made with natural rubber latex.

Indications:

The EASYGRIP FLO-41 System is indicated for delivering FLOSEAL Hemostatic Matrix to bleeding sites through a 5mm or larger trocar.

Intended Purpose / Intended Use:

The EASYGRIP FLO-41 System is intended for delivering FLOSEAL Hemostatic Matrix to bleeding sites through a 5mm or larger trocar.

Contraindications:

None

Warnings:

- **Single use only. Do not reuse. Do not resterilize.** Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.
- Do not use the EASYGRIP FLO-41 System if the packaging is damaged or open, or if the applicator is damaged, as this could lead to compromised device function.
- **Do not use air as a propellant to extrude hemostatic agent from the EASYGRIP FLO-41 System.**
- Do not use the EASYGRIP FLO-41 System to manipulate or retract organs or tissue.
- Do not bend the malleable tip against organs or tissue.

Precautions:

Read the hemostatic agent's instructions prior to using this device.

How Supplied:

The EASYGRIP FLO-41 System is supplied sterile in a sealed plastic tray packed in a product box. The tray contains one (1) applicator device and one (1) empty 1.5 mL syringe. One (1) shipping carton contains six (6) EASYGRIP FLO-41 Systems


NOTE: The provided syringe is used for optional residual product application; do not use prior to complete application of hemostatic agent.

Directions for Use:


Assembly:

1. Prepare the hemostatic agent per its instructions for use.
2. Inspect the integrity of the EASYGRIP FLO-41 System packaging. If the tray, tray seal, or contents have been opened or damaged, do not use.
3. Using aseptic technique, open tray lid where indicated and transfer contents into the sterile field.
4. Attach the hemostatic agent syringe to the EASYGRIP FLO-41 System fill port, taking care to not overtighten.
5. Taking care not to pull the device trigger, fill the reservoir completely with hemostatic agent.

NOTE: Syringe may be removed if desired.

6. Prime cannula by repeatedly pulling trigger (approximately 2-3 pulls) until hemostatic agent is visible in clear malleable tip. The EASYGRIP FLO-41 System is now ready to use.


NOTE: Once primed, take care not to depress trigger prior to application.

7. Ensure the malleable tip is straight before introducing the cannula into a trocar.

Initial Product Application:

8. Carefully approximate the distal end of the malleable tip to the treatment site; use appropriate endoscopic instruments to bend the malleable tip as required for application.
 9. Pull the trigger of the EASYGRIP FLO-41 System to apply the hemostatic agent per its instructions. Multiple pulls may be required to apply desired amount of hemostatic agent.
- NOTE:** For ease of application, apply the hemostatic agent within two (2) hours of reconstitution of the hemostatic agent.

Repeat Product Application (as needed):


1. If additional hemostatic agent is required during the same surgical procedure:
 - load remaining hemostatic agent, taking care not to depress trigger, and fill the reservoir with additional product
 - or-
 - attach new syringe of hemostatic agent

Priming is not needed.

2. Repeat steps 8 and 9 for product application.
- NOTE:** Do not connect the empty syringe or use the saline prior to complete application of the hemostatic agent.

Residual Product Application (optional):

1. If application of residual hemostatic agent is needed, fill the provided 1.5 mL syringe with **non-heparinized saline**.
2. Attach saline syringe to EASYGRIP FLO-41 System and fill the reservoir with only 1.5 mL of non-heparinized saline.

3. Pull the trigger of the EASYGRIP FLO-41 System to apply residual hemostatic agent.
4. If, after expelling all residual product from cannula, an additional hemostatic agent application is required, the cannula must be removed from the patient and purged of saline.

Straighten the malleable tip using appropriate endoscopic instruments and remove the device from the trocar.

5. Perform steps 4-6 from the Assembly section.
- NOTE:** Priming after a saline purge will initially cause saline to be expelled from the cannula; continued priming will fill the cannula with hemostatic agent.

For Disposal:

1. After application of the hemostatic agent straighten the malleable tip with appropriate endoscopic instruments prior to removing the applicator from the trocar.
2. Dispose of the EASYGRIP FLO-41 System according to local regulations. **Do not resterilize.**

Storage:

Store between 2°C (36°F) - 25°C (77°F). Store in original packaging.

Symbols Used in Labeling:

Symbol	Symbol Title and Description	Symbol Reference
	Consult instructions for use; Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2016, 5.4.3
	Do not use if package is damaged; Indicates a medical device that should not be used if the package has been damaged or opened.	ISO 15223-1:2016, 5.2.8
	Sterilized using irradiation; Indicates a medical device that has been sterilized using irradiation.	ISO 15223-1:2016, 5.2.3
	Do not re-sterilize; Indicates a medical device that is not to be resterilized.	ISO 15223-1:2016, 5.2.6
	Temperature limit; Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1:2016, 5.3.7
	Do not re-use; Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1:2016, 5.4.2
	Not made with natural rubber latex; Indicates no presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device.	ISO 15223-1:2016, 5.4.5 and Annex B, Section B2
	Catalogue number; Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1:2016, 5.1.6
	Legal Manufacturer; Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC.	ISO 15223-1:2016, 5.1.1
	Batch code; Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2016, 5.1.5
	Use-by date; Indicates the date after which the medical device is not to be used.	ISO 15223-1:2016, 5.1.4
	Fragile, handle with care; Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1:2016, 5.3.1
	Keep away from sunlight; Indicates a medical device that needs protection from light sources.	ISO 15223-1:2016, 5.3.2
	Keep dry; Indicates a medical device that needs to be protected from moisture.	ISO 15223-1:2016, 5.3.4
	This way up; To indicate correct upright position of the transport package.	ISO 7000 - 0623

U.S. Only: Rx Only

U.S. Only: 1-888-229-0001

Other Countries: Contact your local Baxter representative.

An affiliate of Baxter Healthcare S.A.:

Baxter Healthcare Corporation
One Baxter Parkway
Deerfield, IL 60015, USA

Legal Manufacturer

Baxter Healthcare S.A.
8010 Zurich
Switzerland

0123

Baxter, EASYGRIP FLO-41, FLO-41, and FLOSEAL are trademarks of Baxter International Inc.

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Přidružený subjekt společnosti Baxter Healthcare S.A.:
Baxter, EASYGRIP FLO-41 a FLO-41, FLOSEAL jsou ochranné známky společnosti Baxter International Inc.