

Table of Contents

(continued)

Infusion Systems

Colleague® Volumetric Infusion Devices

Colleague® CX Volumetric Infusion Device
with Guardian Feature150

Colleague® MonoChrome Volumetric
Infusion Device151

Flo-Gard® Volumetric Infusion Devices

Flo-Gard® 6000 Series Volumetric
Infusion Devices.....152

AITECS 2015 Syringe Infusion System

AITECS 2015 Syringe Infusion System152

Table of Contents

(continued)

Infusion Systems

Syringe Pump Sets

Micro-Volume® Extension Sets	153
High Flow Rate Extension Sets	154
Auto Syringe® Sets (Sub-Q-Sets)	154

I.V. Equipment Stands

I.V. Equipment Stands	155
-----------------------------	-----

InfusO.R.™ Pump

InfusO.R.™ Pump	156
InfusO.R.™ SMART™ Labels.....	157
InfusO.R.™ Pump Accessories.....	158

InfusO.R.™ and AutoSyringe® Pump Sets

InfusO.R.™ Pump Sets	159
Economy Sets	159

INFUSION SYSTEMS

Advancing the best in healthcare....for you.

Baxter offers comprehensive clinical inservice support and biomedical training, a nationwide network of technical support and flexible acquisition and financing programs to ensure that every requirement is met. In a changing healthcare environment, Baxter remains your source of support.

For product assistance, please call our Product Information Centre at 1.800.933.0303. This centre is staffed Monday through Friday, 8:00 a.m. to 5:00 p.m. Central Time. If you leave an after-hours message, one of our knowledgeable associates will be happy to return your call on the next business day.

Colleague® Colleague® CX Volumetric Infusion System with Guardian Feature

catalogue number	description	units / case
2M8161	Colleague® CX Volumetric Infusion Pump Single Channel, Colour Screen	1
2M8163	Colleague® CX3 Volumetric Infusion Pump Triple Channel, Colour Screen	1

Features Include:

Guardian Feature & Label Library	Select from 64 drug labels and 32 additional user labels
Uses Standard Set Tubing	Eliminates the need for expensive specialty administration sets
Volume History	Provides detailed volume delivery history
Battery Charge Level	Battery charge level indicator displays battery charge state
Auto Lock / Panel Lock	Locks out unauthorized I.V. set removal or programming
Program Memory	Retains infusion settings for 5 hours after pump is turned off
Flow Rate Range on Each Pump Channel	Primary 0.1 – 1,200 mL / hr Piggyback 0.1 – 500 mL / hr
Volume to be Infused	0.1 mL to 9,999 mL
KVO Flow Rate	0.1 mL increments
Computer Monitoring	
Configurable Downstream Occlusion Values	



catalogue number	description	units / case
2M8151	Colleague® Volumetric Infusion Pump Single-Channel all-purpose infusion device that uses Baxter standard I.V. administration sets.	1
DNM8151	Colleague® Volumetric Infusion Pump Single-Channel - FRENCH	1
2M8153	Colleague® Volumetric Infusion Pump Triple-Channel all-purpose infusion device that uses Baxter standard I.V. administration sets.	1
DNM8153	Colleague® Volumetric Infusion Pump Triple-Channel – FRENCH	1

Features Include:

Micro and Macro Flow Rates
Basic Programming
Advanced Dose Calculation Modes
Piggybacking of Secondary Medications
Automatic Tube Loading
Label Library
Programmable Personalities™ - Storage of up to eight user-defined parameter sets



catalogue number	description	units / case
2D0300	Triple-channel all-purpose Colleague® Infusion Device Syringe Holder (compatible with B.D., Sherwood and Terumo Syringes)	20
2C0169	Colleague Pump Syringe Adapter Set Injection site, non-DEHP tubing, Luer lock adapter Length = 92" (2.3 m) Volume = 4.5 mL	48



non-DEHP

Flo-Gard®

Flo-Gard® 6000 Series Volumetric Infusion System


catalogue number	description	units / case
2M8063	Flo-Gard® 6201 Volumetric Infusion Pump Single-Channel all-purpose infusion device that uses Baxter standard I.V. administration sets.	1
		
2M8064	Flo-Gard® 6301 Volumetric Infusion Pump Dual-Channel all-purpose infusion device that uses Baxter standard I.V. administration sets.	1
		

Features Include:

Automatic Piggybacking
Upstream and Downstream Occlusion Detectors
Air-in-line Detection
Automatic Re-start after an Occlusion Clears
Hardware Based and I.V. Set Free-Flow Protection





AITECS

AITECS 2015 Syringe Infusion System

catalogue number	description	units / case
FAAM3047CA	AITECS 2015 Syringe Infusion System Single Channel Syringe Infusion System that uses Baxter standard syringe administration sets.	1
		

Syringe Pump Sets

Micro-Volume Extension Sets

catalogue number	description	approx overall length	approx volume	units / case
2N3345	MICRO-VOLUME Extension Set with Luer lock adapter Sterile peel pouch	36" (91.4 cm)	0.33 mL	50
		(non-DEHP)		
2N3348	MICRO-VOLUME Extension Set with Luer lock adapter Sterile peel pouch	60" (1.5 m)	0.50 mL	50
		(non-DEHP)		
2N3350	MICRO-VOLUME Extension Set with 0.22 micron air-eliminating, high pressure, pediatric filter, Luer lock adapter Sterile peel pouch	60" (1.5 m)	0.88 mL	50
		(non-DEHP)		
2N3351	MICRO-VOLUME Extension Set with Luer lock adapter, INTERLINK Lever lock cannula Sterile peel pouch	60" (1.5 m)	0.62 mL	50
		(non-DEHP) Interlink IV ACCESS SYSTEM		

Syringe Pump Sets

High Flow Rate Extension Set

catalogue number	description	approx overall length	approx volume	units / case
2N3349	High Flow Rate Extension Set with Luer lock adapter Sterile peel pouch	60" (1.5 m)	2.1 mL	50



non-DEHP

Auto Syringe® Sets (Sub-Q-Sets)

catalogue number	description	approx overall length	approx volume	units / case
1M8450	MICRO-VOLUME Infusion Set with 27 gauge x 1.5 cm (.6") needle Sterile peel pouch	25" (63.5 cm)	0.2 mL	30



1M8451	MICRO-VOLUME Infusion Set with 27 gauge x 1.5 cm (.6") needle Sterile peel pouch	43" (1.1 m)	0.32 mL	30
--------	--	-------------	---------	----



1M8472	SUB-Q-SET Subcutaneous Infusion Set with 27 gauge x 1.2 cm (.475") needle Sterile peel pouch	24" (61.0 cm)	0.23 mL	30
--------	--	---------------	---------	----






1M8463	SUB-Q-SET Subcutaneous Infusion Set with 27 gauge x 1.2 cm (.475") needle Sterile peel pouch	42" (1.1 m)	0.32 mL	30
--------	--	-------------	---------	----



I.V. Stands


I.V. Equipment Stands

catalogue number	description		units / case
2M8953	Five-legged, polished aluminium base, 25" (61 cm) diameter with 2.5" (6 cm) twin wheel casters, Weight 17 lbs. (8 kg)		1
2M8038	Five-legged, polished aluminium base, 25" (61 cm) diameter with 2" (5.1 cm) hooded casters, patient support wheel, Weight 16 lbs. (7 kg)		1
2M8036	Six-legged, heavy base with white epoxy finish, 3" (7 cm) ball bearing casters, 4-outlet UL listed electrical strip, Weight 50 lbs. (23.5 kg)		1

InfusO.R.™ Pump

InfusO.R.™ Pump

catalogue number	description	units / case
2L3100	<p>The InfusO.R.™ Pump is a multi-drug syringe pump for the infusion of anaesthetics, muscle relaxants, analgesics and cardiovascular support drugs in the operating room. The pump is easily programmed by attaching the appropriate drug-specific SMART™ Label to the face of the pump.</p> <p>Pump comes with: 1 SMART™ Label (your choice) Operator's Manual Quick Reference Guide</p>	1



**InfusO.R.™
SMART™
Labels**

InfusO.R.™ Pump SMART™ Labels

catalogue number	description	units / case
2L3226	Alfentanil (Alfenta)	1
2L3246	Amrinone (Inocor)	1
2L3228	Atracurium (Tracrium)	1
2L3239	Atracurium , Peadiatric Dose	1
2L3260	Cisatracurium	1
2L3229	Dobutamine (Dobutrex)	1
2L3230	DOPamine (Intropin)	1
2L3231	Epinephrine (Adrenalin)	1
2L3233	Esmolol (Brevibloc)	1
2L3243	Fentanyl (Sublimaze)	1
2L3241	General mL / hr	1
2L3225	General mL / hr - Bolus	1
2L3232	Isoproterenol (Isuprel)	1
2L3247	Lidocaine (Xylocaine)	1
2L3244	Midazolam (Versed)	1
2L3248	Mivacurium (Mivacron)	1
2L3236	Nitroglycerin (Nitrostat)	1
2L3234	Nitroprusside (Nipride)	1
2L3251	Norepinephrine (Levophed)	1
2L3245	Phenylephrine (Neo-Synephrine)	1
2L3242	Propofol (Diprivan)	1
2L3259	Propofol, Paediatric Dose	1

**InfusO.R.TM
SMARTTM
Labels**


InfusO.R.TM Pump SMARTTM Labels

catalogue number	description	units / case
2L3268	Remifentanil (Ultiva TM) M.A.C. (Monitored Anesthesia Care)	1
2L3269	Remifentanil (Ultiva TM) General Anesthesia	1
2L3237	Succinylcholine (Anectine)	1
2L3235	Sufentanil 50 mcg / MI (Sufenta)	1
2L3249	Sufentanil 5 mcg / mL (Sufenta)	1
2L3240	Sufentanil, Paediatric Dose	1
2L3250	Thopental (Pentothal)	1
2L3227	Vecuronium (Norcuron)	1
2L3238	Vecuronium, Paediatric Dose	1

**InfusO.R.TM
Pump
Accessories**




InfusO.R.TM Pump Accessories

catalogue number	description	units / case
2L3223	Sleeve to hold SMART TM Labels	1
2L3212	Pole Clamp, Non-Locking	1



**InfusO.R.™
and Auto
Syringe®
Pump Sets**

InfusO.R.™ Pump Sets

catalogue number	description	approx overall length	approx priming volume	approx inner diameter	units / case
2C9215	Anti-siphon Anesthesia Set with BLUE stripe for easy identification, anti-siphon valve prevents inadvertent drug delivery during syringe changes Sterile peel pouch	72" (1.8 m)	1.6 mL	0.040"	60
					
2C9216	Anesthesia Set with GREEN stripe for easy identification, Luer lock adapter	72" (1.8 m)	1.6 mL	0.040"	60
					
2C9218	Universal Set with YELLOW stripe for easy identification, larger bore for use with drugs of higher viscosity, Luer lock adapter	72" (1.8 m)	2.8 mL	0.055"	60
					

Economy Sets

catalogue number	description	approx overall length	approx priming volume	approx inner diameter	units / case
2C9219	Extension Set, General All-Purpose Set	60" (1.5 m)	1.4 mL	0.040"	60
2C9221	Extension Set, High Flow Rate larger bore for use with drugs of higher viscosity	60" (1.5 m)	2.5 mL	0.055"	60