

Table of Contents

(continued)

Pain Management

Syndeo PCA Syringe Pump.....	162
Syndeo Pump Sets	162
Ipump™ Pain Management System	164
IPump™ Sets for use with IPump™ Bags....	164
IPump™ Sets for use with I.V. Bags.....	165
IPump™ Bags	167
IPump™ Accessories	167

PAIN MANAGEMENT

Advancing the best in healthcare....for you.

Baxter provides PCA / PCEA Pain Management Systems that streamline pain therapy with a combination of intelligent design and simple operation.

Our goal is to consistently provide the highest quality pain management products, technical information, clinical expertise and consulting services.


Sydeo PCA Syringe Pump

Sydeo PCA Syringe Pump

catalogue number	description	units / case
2L3113	SYNDEO PCA Syringe Pump	1

Pain management can be complicated but that doesn't mean your PCA infusion system should be. The new SYNDEO PCA Syringe Pump streamlines pain therapy with a powerful combination of intelligent design and simple operations.

SYNDEO



Sydeo Pump Sets





Sydeo Pump Sets

catalogue number	description	approx overall length	approx priming volume	approx inner diameter	units / case
2L3506	8" Anti-Reflux Y-Set Microbore tubing, removable non-latex injection port for intermittent injections or connection of additional line with backcheck valve	8" (20 cm)	1.4 mL	0.06"	50
2C9205	Extension Set Tamper resistant, anti-siphon PCA	60" (1.5 m)	1.2 mL	0.03"	50



**Sydeo
Pump Sets**


Sydeo Pump Sets

catalogue number	description	approx overall length	approx priming volume	approx inner diameter	units / case
2C9206	Extension Set Tamper resistant, anti-siphon PCA	96" (2.4 m)	1.6 mL	0.03"	50
					
2L3508	Extension Set Anti-siphon "Y" with anti-reflux	98" (2.5 m)	2.4 mL	0.03"	50
					
2L3507	Combination Set, Tamper resistant, anti-siphon PCA with anti-reflux Y-set	101" (2.6 m)	2.5 mL	0.03"	50
					
2L3509	Extension Set Anti-siphon "Y" with anti-reflux	102" (2.6 m)	3.3 mL	0.03"	50
					

IPump™ Pain Management System



IPump™ Pain Management System

catalogue number	description	units / case
2L3107	The IPump™ Pain Management System answers the day-to-day needs of the clinical environment with one-on-one flexibility. From programming to data retrieval, the IPump™ System puts customizable pain management right where you want it.	1



IPump™ Sets

IPump™ Sets for use with IPump™ Bags

catalogue number	description	approx overall length	approx priming volume	tubing inner diameter	units / case
2L3510	Anti-Siphon Luer Lock Set Attaches to medication bag with Luer lock connector, anti-siphon valve minimizes potential for free-flow.	73" (1.8 m)	1.5 mL	0.030"	20
					
2L3511	Anti-Siphon Luer Lock Set Attaches to medication bag with Luer lock connector, anti-siphon valve minimizes potential for free-flow.	102" (2.5 m)	1.8 mL	0.030"	20
					
2L3504	EPIDURAL Anti-Siphon Luer Lock Set Attaches to medication bag with Luer lock connector, anti-siphon valve minimizes potential for free flow, YELLOW stripe denotes EPIDURAL use.	104" (2.6 m)	1.6 mL	0.030"	20

IPump™ Sets for use with IPump™ Bags



catalogue number	description	approx overall length	approx priming volume	tubing inner diameter	units / case
2L3520	Air-Eliminating Anti-Siphon Luer Lock Set Attaches to medication bag with Luer lock connector, air-eliminating 0.22 micron filter , anti-siphon valve minimizes potential for free-flow.	106" (2.7 m)	2.1 mL	0.030"	20
2L3526	Air-Eliminating Anti-Siphon Luer Lock Set with Integral Y-Site Attaches to medication bag with Luer lock connector, air-eliminating 0.22 micron filter , anti-siphon valve minimizes potential for free-flow, integral anti-reflux Y-site connection.	106" (2.7 m)	3.2 mL	0.030"	20

IPump™ Sets for use with I.V. Bags

catalogue number	description	approx overall length	approx priming volume	tubing inner diameter	units / case
2L3506	8" Anti-Reflux Y-Set Microbore tubing, removable non-latex injection port for intermittent injections or connection of additional line with backcheck valve.	8" (20 cm)	1.4 mL	0.061"	50
2L3525	Air-Eliminating Anti-Siphon Spike Set with Integral Y-Site Attaches to medication bag with spike port, air-eliminating 0.22 micron filter , anti-siphon valve minimizes potential for free-flow, integral anti-reflux Y-site connection.	106" (2.7 m)	4.1 mL	0.030"	20




IPump™ Sets for use with I.V. Bags

catalogue number	description	approx overall length	approx priming volume	tubing inner diameter	units / case
2L3512	EPIDURAL Anti-Siphon Spike Set Attaches to medication bag with spike port, anti-siphon valve minimizes potential for free-flow, YELLOW stripe denotes EPIDURAL use.	107" (2.7 m)	2.4 mL	0.030"	20
					
2L3513	Air-Eliminating Anti-Siphon Spike Set Attaches to medication bag with spike port, air-eliminating 0.22 micron filter , anti-siphon valve minimizes potential for free-flow.	107" (2.7 m)	2.4 mL	0.030"	20
					
2L3523	Air-Eliminating Anti-Siphon Filter Spike Set for 500 mL Bag Cover Attaches to medication bag with administration port, 16" upper segment to accommodate the 500 mL bag cover, air-eliminating 0.22 micron filter .	120" (3.0 m)	2.5 mL	0.030"	20
2L3522	EPIDURAL Anti-Siphon Spike Set for use with 500 mL Bag Cover Attaches to medication bag with administration port, 16" upper segment to accommodate the 500 mL bag cover, YELLOW stripe on the lower segment denotes EPIDURAL use.	121" (3.1 m)	2.5 mL	0.030"	20
2L3527	Air-Eliminating Anti-Siphon Integral Spike with Y-Site For use with 500mL bag cover	126" (3.2 m)	4.2 mL	0.030"	20

IPump™ Bags

IPump™ Bags




catalogue number	description	units / case
2L3256	100 mL Drug Bag Locks inside 100 mL bag cover Female Luer lock adapter	10
2L3257	250 mL Drug Bag Locks inside 250 mL bag cover Female Luer lock adapter	10



IPump™ Accessories






IPump™ Accessories

catalogue number	description	units / case
2L3400	Printer Adapter Adapts IPump™ to standard printers (requires printer adapter cable)	1
2L3211	Locking Pole Clamp Includes 2 keys	1
2L3212	Non-Locking Pole Clamp	1
2L3112	Configuration Cable Allows pump configuration to be transferred electronically from one IPump™ to another IPump™.	1



**IPump™
Accessories**

IPump™ Accessories

catalogue number	description	units / case
2L3210	AC Adapter 	1
2L3218	100 mL Bag Cover	1
2L3220	250 mL Bag Holder	1
2L3217	250E Bag Cover 	1
2L3221	500 mL Bag Cover 	1
2L3219	Carrying Case 	1
2L3402	Printer Adapter Cable Allows connection of Printer Adapter to IPump™ printer port. 	1