

Peritoneal Dialysis

PD Solutions Designed to optimize patient prescription



Dianeal®

For acute or chronic renal failure patients

PD4: Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

PD101: Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

Available in 0.5, 1.5, 2.5 and 4.25% dextrose



Nutrineal® PD4

For treatment of Protein Malnutrition in Peritoneal Dialysis Patients, in a Glucose-Free, Phosphate-Free Exchange.

Na 132, Ca 1.25, Mg 0.25, Cl 105, Lactate 40 mmol/L



Extraneal®

Solution for the Long Dwell in CAPD and APD Patients, using Icodextrin instead of Glucose as the Primary Osmotic Agent.

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L



Physioneal® 40

A Solution that provides a Physiologically Balanced Environment for peritoneal membrane.

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15mmol/L, Bicarb 25mmol/L

Available in 1.36, 2.27 and 3.86% glucose

Dianeal® PD4, PD101 and Physioneal® 40 Prescribing Information available upon request.

Nutrineal® PD4 (1.1% Amino Acid Peritoneal Dialysis Solution) and Extraneal® (7.5% Icodextrin) Product Monographs available upon request.

Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

TWIN BAG®



The system of choice for Canadian CAPD patients.

When it comes to your patients, trust TwinBag

- Simple Design, Successful Therapy
- Engineered for optimal flush efficiency
- Protects against touch contamination
- Allows effective removal of potential contaminants
- Designed to minimize fibrin blockage
- Latex free

CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	DRAIN BAG VOLUME (mL)	PACK FACTOR
DIANEAL® - PD4				
Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L				
JB9415	1.50% Dextrose	1500	2000	6/case
JB9425	2.50% Dextrose	1500	2000	6/case
JB9435	4.25% Dextrose	1500	2000	6/case
JB9416	1.50% Dextrose	2000	2000	6/case
JB9426	2.50% Dextrose	2000	2000	6/case
JB9436	4.25% Dextrose	2000	2000	6/case
JB9418	1.50% Dextrose	2500	3000	4/case
JB9428	2.50% Dextrose	2500	3000	4/case
JB9438	4.25% Dextrose	2500	3000	4/case
JB9413	1.50% Dextrose	3000	3000	4/case
JB9423	2.50% Dextrose	3000	3000	4/case


DIANEAL® - PD101				
Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L				
JB9806	0.50% Dextrose	2000	2000	6/case
JB9816	1.50% Dextrose	2000	2000	6/case
JB9826	2.50% Dextrose	2000	2000	6/case
JB9836	4.25% Dextrose	2000	2000	6/case
JB9818	1.50% Dextrose	2500	3000	4/case
JB9828	2.50% Dextrose	2500	3000	4/case

Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

TWIN BAG®

(Continued)



COLOUR CODED CAP CONFIRMS DEXTROSE CONCENTRATION

- White = 0.50% dextrose (Dianeal®)
- Yellow = 1.50% dextrose (Dianeal®) and 1.36% glucose (Physioneal®)
- Green = 2.50% dextrose (Dianeal®) and 2.27% glucose (Physioneal®)
- Orange = 4.25% dextrose (Dianeal®) and 3.86% glucose (Physioneal®)
- Blue = 1.1% amino acids (Nutrineal®)
- Purple = 7.5% Icodextrin (Extraneal®)

CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	DRAIN BAG VOLUME (mL)	PACK FACTOR
------------------	---------------	------------------	-----------------------	-------------



NUTRINEAL®

Na 132, Ca 1.25, Mg 0.25, Cl 105, Lactate 40 mmol/L

JB2019	1.1% Amino Acid Solution	2000	2000	6/case
JB2005	1.1% Amino Acid Solution	2500	3000	4/case



EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

JB9912	7.5% Icodextrin Solution	2000	2000	6/case
JB9913	7.5% Icodextrin Solution	2500	3000	4/case



PHYSIONEAL®

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25mmol/L

JB8249	1.36% Glucose	2000	2500	6/case
JB8251	2.27% Glucose	2000	2500	6/case
JB8254	3.86% Glucose	2000	2500	6/case
JB8250	1.36% Glucose	2500	3000	5/case
JB8252	2.27% Glucose	2500	3000	5/case

Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
------------------	-------------	-------------

Disconnect Y-Sets

5C4366P	CAPD Disconnect Y-Set, Bulk	60/case
---------	-----------------------------	---------

Miscellaneous CAPD Supplies

5C4462	Empty Sterile Bag, 3L, Spike	48/case
XMC4284	Empty Sterile Heater/Drain Bag, 3L, Luer Lock	30/box

CAPD Solution Warmer



CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
JBWarmerU	CAPD Solution Warmer	1 each

Designed to provide patients with simple and efficient system for warming TwinBag solutions.

- Ensures integrity of PD solutions
- Preset to a comfortable temperature
- Secondary overheat thermal protector shutoff
- Lightweight and portable
- Requires only 20 watts of electricity
- Made in Canada



Peritoneal Dialysis

AMIA Automated Peritoneal Dialysis System with Sharesource Connectivity Platform

AVAILABLE FOR ADULT PATIENTS

- 40% smaller and 30% lighter than the HOMECHOICE PRO Automated PD System
- Provides a user-friendly PD experience that guides patients through a treatment or therapy
- Equipped with voice guidance, touch screen, animations, and remote monitoring capabilities



DEVICES

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
5C9320	AMIA APD SYSTEM V2.0	1
5C6012	Country Kit, V 2.0, Canada, English	1
5C9346	Cycler Travel Case	1
5C9400	AMIA Sharesource Platform	NA
R5C5522	Personal Data Key 5-Pack	5
5C6M20	Kit cellular Modem and Power	1
5C5479	Set 4 prong. APD Cycler	30

Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

HomeChoice Pro™ and HomeChoice™ Automated Peritoneal Dialysis Systems


**AVAILABLE FOR
PEDIATRIC TO ADULT
PATIENTS**



- Pneumatic pumps – no scales!
- No re-calibration necessary after moving the machine
- Patient management with ProCard and RenalSoft PD

The low-fill mode option for pediatric and low-fill volume adult patients offers:

- Fill volume range from 60 – 3000 mL
- Programmable minimum drain time and volume parameters for more closely regulated therapy
- Alarm limits and programmable positive and negative UF alarms ensure proper drainage and increase treatment flexibility



DEVICES

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
5C4471R	HomeChoice™ Automated PD System	1 each
5C8310P	HomeChoice PRO™ Automated PD System (110 volt)	1 each

Peritoneal Dialysis Software



- APD therapies are written by the clinician and transferred to the cyclor with the PRO Card
- Window's based software enables increased monitoring of patient adherence, drain function and ultrafiltration

R5C4522	HomeChoice PRO™ Card	6/case
---------	----------------------	--------

Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

APD Luer Lock Solutions

CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	PACK FACTOR
------------------	---------------	------------------	-------------

DIANEAL® - PD4

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

JB4769L	1.50% Dextrose	3000	4/case
JB4770L	2.50% Dextrose	3000	4/case
JB4771L	4.25% Dextrose	3000	4/case
JB9719L	1.50% Dextrose	5000	2/case
JB9729L	2.50% Dextrose	5000	2/case
JB9749L	4.25% Dextrose	5000	2/case

DIANEAL® - PD101

Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

JB4821L	0.50% Dextrose	2000	6/case
JB5618L	1.50% Dextrose	2500	5/case
JB5658L	2.50% Dextrose	2500	5/case
JB5619L	1.50% Dextrose	5000	2/case
JB5659L	2.50% Dextrose	5000	2/case
JB5669L	4.25% Dextrose	5000	2/case

NUTRINEAL®

Na 132, Ca 1.25, Mg 0.25, Cl 105, Lactate 40 mmol/L

JB2006LP	1.1% Amino Acid Solution	2500	5/case
----------	--------------------------	------	--------

EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

JB9923LP	7.5% Icodextrin Solution	2500	5/case
----------	--------------------------	------	--------

PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L

JB8291	1.36% Glucose	2500	5/case
JB8290	2.27% Glucose	2500	5/case
JB8289	3.86% Glucose	2500	5/case

PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L

CJPE8280	1.36% Glucose	5000	2/case
CJPE8281	2.27% Glucose	5000	2/case
CJPE8282	3.86% Glucose	5000	2/case

Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

APD DISPOSABLE LUER SETS

L5C4531	3-Prong Cyclor Set with Cassette	30/case
R5C4479C	4-Prong Cyclor Set with Cassette	30/case
R5C4478	8-Prong Cyclor Set with Cassette	30/case
R5C4427	5-Prong Manifold Set	30/case
N5C8305C	Low Recirculation Volume APD Set	30/case

- Low Recirculation Volume APD Set has in-line medication port, LuerLock to bag connectors, a shrouded patient connector which covers all transfer set threads, and a smaller internal diameter and shorter patient line length which reduces patient line recirculation volume to 17mL for a more efficient therapy.

APD EXTENSION LINE



5C4480C	12 Foot Patient Extension Line, Easy Lock™ Connector (3.6M)	30/case
---------	---	---------

APD DRAIN SUPPLIES

XMC4284	Empty Sterile Heater\Drain Bag, 3L, Luer Lock	30/box
5C4462	Empty Sterile Bag, 3L, Spike	48/case
5C4476	Effluent Sample Bag, 150ml, Luer	6/case
5C4145P	Cyclor Drainage Bag Set (2 x 15L per pack), Spike	15 pk/case
5C4512	APD Drain Manifold (Connects 2 15 Litre drain bags), Spike	30/case
5C4464P	12 Foot Drain Line Extension Set (3.6M), Spike	30/case

APD CARRYING CASES

5C4406	HomeChoice: Soft-sided Carrying Case	1 each
5C8346	HomeChoice: Shipping Case with wheels	1 each
5C4971	Quiet Coil (ON HOLD)	1 each



Peritoneal Dialysis

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
Transfer Sets		
Disconnect System		
5C4483	MiniCap Extended Life Transfer Set with Twist Clamp (10cm)*	6/case
5C4482	MiniCap Extended Life Transfer Set with Twist Clamp (15cm)*	6/case
R5C4484	MiniCap Extended Life Transfer Set with Twist Clamp (22cm)*	6/case
5C4449	MiniCap Extended Life Transfer Set with Roller Clamp (15cm)*	6/case
Non-Disconnect System		
5C4160	System III Transfer Set with Spike Connector (120cm)*	12/case
* Extended Life Transfer Set Wear – 6 months		
Transfer Set Adapters		
5C4136	Luer-Lock to Spike Adapter Set	12/case
Catheter Supplies		
5C4129	Locking Titanium Adapter (2 pieces)	5/bx,24bx/ca
5C4169	Locking Cap for 5C4129 Titanium Adapter	25/case
PC4171	Catheter Clamp	12/case
5C4319	Immobile™ catheter and transfer set stabilizing device	30/case
5C4320	Immobile A/C, Catheter stabilizer	30/case
5C4381	Immobile™ catheter stabilizing device	30/case
Disconnect Caps with Pre-filled Povidone-Iodine Solution		
5C4466P	MiniCap Disconnect Cap (English labeling only)	60/case
SPC4466	MiniCap Disconnect Cap (same as RPC4466)	60/case
SPC4486	OptiCap Disconnect Cap	30/case
Clamps		
5C4527	Outlet Port Clamp, Short-Nosed, Blue	12/case
5C4957	UltraClamp, Red (Not for use as outlet port clamp)	12/case
EZ-Aide™		
5C4505	Assist Device for CAPD or APD patients	1 each

Bolded codes are preferred products



5C4483



5C4129



SPC4486



5C4527



5C4957



5C4505