
Peritoneal Dialysis

PD Solutions Designed to optimize patient prescription



Dianeal®

For acute or chronic renal failure patients

PD4: Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

PD101: Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

Available in 0.5, 1.5, 2.5 and 4.25% dextrose



Extraneal®

Solution for the Long Dwell in CAPD and APD Patients, using Icodextrin instead of Glucose as the Primary Osmotic Agent.

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L



Physioneal® 40

A Solution that provides a Physiologically Balanced Environment for peritoneal membrane.

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15mmol/L, Bicarb 25mmol/L

Available in 1.36, 2.27 and 3.86% glucose

Dianeal® PD4, PD101 and Physioneal® 40 Prescribing Information available upon request.
Extraneal® (7.5% Icodextrin) Product Monographs available upon request.

Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

TWIN BAG®

The system of choice for Canadian CAPD patients.

When it comes to your patients, trust TwinBag

- Simple Design, Successful Therapy
- Engineered for optimal flush efficiency
- Protects against touch contamination
- Allows effective removal of potential contaminants
- Designed to minimize fibrin blockage
- Latex free



| CATALOGUE NUMBER | CONCENTRATION | FILL VOLUME (mL) | DRAIN BAG VOLUME (mL) | PACK FACTOR |
|------------------|---------------|------------------|-----------------------|-------------|
|------------------|---------------|------------------|-----------------------|-------------|

DIANEAL® - PD4

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

| | | | | |
|--------|----------------|------|------|--------|
| JB9415 | 1.50% Dextrose | 1500 | 2000 | 6/case |
| JB9425 | 2.50% Dextrose | 1500 | 2000 | 6/case |
| JB9435 | 4.25% Dextrose | 1500 | 2000 | 6/case |
| JB9416 | 1.50% Dextrose | 2000 | 2500 | 6/case |
| JB9426 | 2.50% Dextrose | 2000 | 2500 | 6/case |
| JB9436 | 4.25% Dextrose | 2000 | 2500 | 6/case |
| JB9418 | 1.50% Dextrose | 2500 | 3000 | 4/case |
| JB9428 | 2.50% Dextrose | 2500 | 3000 | 4/case |
| JB9438 | 4.25% Dextrose | 2500 | 3000 | 4/case |

DIANEAL® - PD101

Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L


| | | | | |
|--------|----------------|------|------|--------|
| JB9806 | 0.50% Dextrose | 2000 | 2500 | 6/case |
| JB9816 | 1.50% Dextrose | 2000 | 2500 | 6/case |
| JB9826 | 2.50% Dextrose | 2000 | 2500 | 6/case |
| JB9836 | 4.25% Dextrose | 2000 | 2500 | 6/case |

Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

TWIN BAG®

(Continued)



COLOUR CODED CAP CONFIRMS DEXTROSE CONCENTRATION

- White = 0.50% dextrose (Dianeal®)
- Yellow = 1.50% dextrose (Dianeal®) and 1.36% glucose (Physioneal®)
- Green = 2.50% dextrose (Dianeal®) and 2.27% glucose (Physioneal®)
- Orange = 4.25% dextrose (Dianeal®) and 3.86% glucose (Physioneal®)
- Purple = 7.5% Icodextrin (Extraneal®)

| CATALOGUE NUMBER | CONCENTRATION | FILL VOLUME (mL) | DRAIN BAG VOLUME (mL) | PACK FACTOR |
|------------------|---------------|------------------|-----------------------|-------------|
|------------------|---------------|------------------|-----------------------|-------------|



EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

| | | | | |
|--------|--------------------------|------|------|--------|
| JB9912 | 7.5% Icodextrin Solution | 2000 | 2500 | 6/case |
| JB9913 | 7.5% Icodextrin Solution | 2500 | 3000 | 4/case |






PHYSIONEAL®

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25mmol/L

| | | | | |
|--------|---------------|------|------|--------|
| JB8249 | 1.36% Glucose | 2000 | 2500 | 6/case |
| JB8251 | 2.27% Glucose | 2000 | 2500 | 6/case |
| JB8254 | 3.86% Glucose | 2000 | 2500 | 6/case |
| JB8250 | 1.36% Glucose | 2500 | 3000 | 5/case |
| JB8252 | 2.27% Glucose | 2500 | 3000 | 5/case |


Peritoneal Dialysis

CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|------------------------------------|---|-------------|---|
| Disconnect Y-Sets | | | |
| 5C4366P | CAPD Disconnect Y-Set, Bulk | 60/case |  |
| Miscellaneous CAPD Supplies | | | |
| 5C4462 | Empty Sterile Bag, 3L, Spike | 48/case |  |
| XMC4284 | Empty Sterile Heater/Drain Bag, 3L, Luer Lock | 30/box |  |

Peritoneal Dialysis

CAPD Solution Warmer

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|------------------|----------------------|-------------|---|
| JBWarmerU | CAPD Solution Warmer | 1 each |  |

Designed to provide patients with simple and efficient system for warming TwinBag solutions.

- Ensures integrity of PD solutions
- Preset to a comfortable temperature
- Secondary overheat thermal protector shutoff
- Lightweight and portable
- Requires only 20 watts of electricity
- Made in Canada

Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

HomeChoice Claria

- Initial Drain (I-Drain) Bypass - Maximizing patient comfort and confidence
- Sharesource 2-way connectivity platform for Remote Patient Management (RPM)
- 2-line OLED screen that eliminates alternating messages and helps improve the user experience (1)
- Increased number of languages to 38



HomeChoice Claria

DEVICES

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR |
|------------------|--|-------------|
| 5C6M40 | HOMECHOICE CLARIA APD | 1 each |
| 5C6M40R | HOMECHOICE CLARIA APD REFURB | 1 each |
| 5C7346 | HOMECHOICE CLARIA CARRYINGCASE | 1 each |
| 5CGM01 | CLARIA SHARESOURCE PLATFORM | 1 each |
| 5C6M70 | 4G KIT CELLULAR MODEM AND POWER | 1 each |
| BXU005345 | POWER CORD FOR CLARIA 125V NEMA | 1 each |
| 0719003471A2* | HOMECHOICE CLARIA PATIENT AT HOME GUIDE (PAHG) | 1 each |
| 0719001577C1** | HOMECHOICE CLARIA Clinician guide English | 1 each |
| 0719003513A1*** | HOMECHOICE CLARIA Clinician guide Addendum English | 1 each |



4G Modem



¹⁾ Baxter International Inc. Deerfield, IL

*0719003471A2FRE: PAHG French. **0719001577C1FRE: Clinician guide French. ***0719003513A1FRE: Clinician guide addendum French.

Peritoneal Dialysis

AMIA Automated Peritoneal Dialysis System with Sharesource Connectivity Platform

AVAILABLE FOR ADULT PATIENTS

- 40% smaller and 30% lighter than the HOMECHOICE PRO Automated PD System
- Provides a user-friendly PD experience that guides patients through a treatment or therapy
- Equipped with voice guidance, touch screen, animations, and remote management capabilities



DEVICES

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR |
|------------------|-------------------------------------|-------------|
| 5C9320 | AMIA APD SYSTEM V2.0 | 1 |
| 5C9320R | AMIA APD SYSTEM V2.0 Refurb | 1 |
| 5C6012 | Country Kit, V 2.0, Canada, English | 1 |
| 5C9346 | Cycler Travel Case | 1 |
| 5C9400 | AMIA Sharesource Platform | NA |
| R5C5522 | Personal Data Key 5-Pack | 5 |
| 5C6M70 | 4G Kit cellular Modem and Power | 1 |
| 5C5479 | Set 4 prong. APD Cycler | 30 |



Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

HomeChoice Pro™ and HomeChoice™ Automated Peritoneal Dialysis Systems

**AVAILABLE FOR
PEDIATRIC TO ADULT
PATIENTS**

- Pneumatic pumps – no scales!
- No re-calibration necessary after moving the machine
- Patient management with ProCard and RenalSoft PD

The low-fill mode option for pediatric and low-fill volume adult patients offers:

- Fill volume range from 60 – 3000 mL
- Programmable minimum drain time and volume parameters for more closely regulated therapy
- Alarm limits and programmable positive and negative UF alarms ensure proper drainage and increase treatment flexibility



HomeChoice Pro™

DEVICES

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR |
|------------------|--|-------------|
| 5C4471RU* | HomeChoice™ Automated PD System | 1 each |
| 5C4471RUL | HomeChoice™ Automated PD System with Mandarin Language | 1 each |
| 5C8310RU** | HomeChoice PRO™ Automated PD System (110 volt) | 1 each |

Peritoneal Dialysis Software

- APD therapies are written by the clinician and transferred to the cyclor with the PRO Card
- Window's based software enables increased monitoring of patient adherence, drain function and ultrafiltration

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|------------------|---|-------------|---------|
| R5C4522 | HomeChoice PRO™ Card | 6/case | |
| 4994000449 | Homechoice Pro USB SmartMedia Card Reader | 1 each | |

*5C4471RU is replacing 5C4471R. **5C8310RU is replacing 5C8310P

Peritoneal Dialysis

APD – AUTOMATED PERITONEAL DIALYSIS

APD Luer Lock Solutions

| CATALOGUE NUMBER | CONCENTRATION | FILL VOLUME (mL) | PACK FACTOR |
|------------------|---------------|------------------|-------------|
|------------------|---------------|------------------|-------------|

DIANEAL® - PD4

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

| | | | |
|---------|----------------|------|--------|
| JB4769L | 1.50% Dextrose | 3000 | 4/case |
| JB4770L | 2.50% Dextrose | 3000 | 4/case |
| JB4771L | 4.25% Dextrose | 3000 | 4/case |
| JB9719L | 1.50% Dextrose | 5000 | 2/case |
| JB9729L | 2.50% Dextrose | 5000 | 2/case |
| JB9749L | 4.25% Dextrose | 5000 | 2/case |

DIANEAL® - PD101

Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

| | | | |
|---------|----------------|------|--------|
| JB4821L | 0.50% Dextrose | 2000 | 6/case |
| JB5618L | 1.50% Dextrose | 2500 | 5/case |
| JB5658L | 2.50% Dextrose | 2500 | 5/case |
| JB5619L | 1.50% Dextrose | 5000 | 2/case |
| JB5659L | 2.50% Dextrose | 5000 | 2/case |
| JB5669L | 4.25% Dextrose | 5000 | 2/case |

EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

| | | | |
|----------|--------------------------|------|--------|
| JB9923LP | 7.5% Icodextrin Solution | 2500 | 5/case |
|----------|--------------------------|------|--------|

PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L

| | | | |
|--------|---------------|------|--------|
| JB8291 | 1.36% Glucose | 2500 | 5/case |
| JB8290 | 2.27% Glucose | 2500 | 5/case |
| JB8289 | 3.86% Glucose | 2500 | 5/case |






PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L







| | | | |
|----------|---------------|------|--------|
| CJPE8280 | 1.36% Glucose | 5000 | 2/case |
| CJPE8281 | 2.27% Glucose | 5000 | 2/case |
| CJPE8282 | 3.86% Glucose | 5000 | 2/case |

Peritoneal Dialysis





APD – AUTOMATED PERITONEAL DIALYSIS

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|---------------------------------|--|-------------|---|
| APD DISPOSABLE LUER SETS | | | |
| R5C4479C | 4-Prong Cyclor Set with Cassette | 30/case |  |
| R5C4478 | 8-Prong Cyclor Set with Cassette | 30/case |  |
| R5C4427 | 5-Prong Manifold Set | 30/case |  |
| N5C8305C | Low Recirculation Volume APD Set Low Recirculation Volume APD Set has in-line medication port, Luer Lock to bag connectors, a shrouded patient connector which covers all transfer set threads, and a smaller internal diameter and shorter patient line length which reduces patient line recirculation volume to 17mL for a more efficient therapy. | 30/case |  |
| APD EXTENSION LINE | | | |
| 5C4480C | 12 Foot Patient Extension Line, Easy Lock™ Connector (3.6M) | 30/case |  |









Peritoneal Dialysis

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|---------------------------|---|-------------|---|
| APD DRAIN SUPPLIES | | | |
| XMC4284 | Empty Sterile Heater/Drain Bag, 3L, Luer Lock | 30/box |  |
| 5C4462 | Empty Sterile Bag, 3L, Spike | 48/case |  |
| 5C4476 | Effluent Sample Bag, 150ml, Luer | 6/case |  |
| 5C4145NS | Cycler Drainage Bag Set (2 x 15L per pack), Spike Non-Sterile | 15 pk/case |  |
| 5C4512 | APD Drain Manifold (Connects 2 15 Litre drain bags), Spike | 30/case |  |
| 5C4464P | 12 Foot Drain Line Extension Set (3.6M), Spike | 30/case |  |






Peritoneal Dialysis

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|---------------------------|---------------------------------------|-------------|---|
| APD CARRYING CASES | | | |
| 5C7346 | Homechoice Claria Carrying Case | 1 each |  |
| 5C9346 | AMIA Cyclor Travel Case | 1 each |  |
| 5C4406 | HomeChoice: Soft-sided Carrying Case | 1 each |  |
| 5C8346 | HomeChoice: Shipping Case with wheels | 1 each |  |

Peritoneal Dialysis

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|---|--|----------------|---|
| Transfer Sets -Disconnect System | | | |
| 5C4483 | MiniCap Extended Life Transfer Set with Twist Clamp (10cm)* | 6/case |  |
| 5C4482 | MINICAP EXTENDED LIFE TRANSFER SET WITH TWIST CLAMP (15CM)* | 6/case | |
| R5C4484 | MiniCap Extended Life Transfer Set with Twist Clamp (22cm)* | 6/case | |
| 5C4449 | MiniCap Extended Life Transfer Set with Roller Clamp (15cm)* | 6/case |  |
| -Non-Disconnect System | | | |
| Transfer Set Adapters | | | |
| 5C4136 | Luer-Lock to Spike Adapter Set | 12/case |  |
| Catheter Supplies | | | |
| 5C4129 | Locking Titanium Adapter (2 pieces) | 5/bx,24bx/ca |  |
| 5C4169 | Locking Cap for 5C4129 Titanium Adapter | 25/case |  |
| PC4171 | Catheter Clamp (same as RPC4171) | 12/case |  |
| 5C4319 | IMMOBILE™ CATHETER AND TRANSFER SET STABILIZING DEVICE | 30/case |  |
| 5C4320 | IMMOBILE A/C, CATHETER STABILIZER | 30/case |  |

Peritoneal Dialysis

| CATALOGUE NUMBER | DESCRIPTION | PACK FACTOR | PICTURE |
|---|---|-------------|---|
| Disconnect Caps with Pre-filled Povidone-Iodine Solution | | | |
| 5C4466P | MiniCap Disconnect Cap (English labeling only) | 60/case |  |
| SPC4466 | MiniCap Disconnect Cap (same as RPC4466) | 60/case | |
| SPC4456 | Flexicap Disconnect Cap | 30/case |  |
| Clamps | | | |
| 5C4527 | Outlet Port Clamp, Short-Nosed, Blue | 12/case |  |
| 5C4957 | UltraClamp, Red (Not for use as outlet port clamp) | 12/case |  |
| EZ-Aide™ | | | |
| 5C4505 | Assist Device for CAPD or APD patients | 1 each |  |

**Bolded codes are preferred products*